



Fact Sheet

O2 Cold and Compression Back Wrap

The O2 Pneumatic/Cold Back wrap is designed to achieve the cold AND compression effects needed to rehabilitate an injury or recover quickly from post-operative care. This desired outcome is achieved without costly and cumbersome flow pumps or ice bags that leak.

People enjoy the ease of O2 application and air as necessary. The O2 Back wrap does not require water or electricity, simply place it back inside the refrigerator when finished.

The O2's anatomic design correctly covers the affected areas.

- ◆ Easy to adjust
- ◆ Wearable over clothing
- ◆ Frostbite proof
- ◆ Quick pump top alleviate pain faster
- ◆ Use for sprains, strains, sports injuries, rehab, post surgery, arthritis, fatigues, muscle tensions, etc.
- ◆ Drip free fabric